



Remote Access Guide for PROXIMA 2A

Version 06 April 2023

Revised 26 May 2023

Outline

- Preparing for your Remote Access Session
 - Sending Dewars
 - Activation of your account for Remote Access
 - Testing your internet connection
- Updates and changes
 - New remote access server
 - New MXCuBE workstation
 - New NoMachine version
 - ADXV & ALBULA
- Connection via NoMachine
 - First steps
 - Custom Display
 - Connection Errors
- Connected
 - PX2-A Overview
 - X-terminals
 - Firefox & Webcams
 - MXCuBE
 - ALBULA & ADXV
- Other stuff
 - HDF5 data architecture
 - Data Processing
 - Process1
 - XDSME
 - Help, Hints, FAQs & Bugs

Preparing for REMOTE ACCESS : Sending Dewars


- Please read the tutorial on the PX2-A web site
 - <https://www.synchrotron-soleil.fr/fr/lignes-de-lumiere/proxima-2a>
- A week before your session, send your Dewar to SOLEIL
 - 2 options:
 - Courier-delivered (Fedex, DHL, UPS,...)
 - Pre-paid Air Way Bills (SOLEIL, iNext or User)
 - Hand-delivered directly to the SOLEIL gate
 - Please organise a single delivery with the local contact
 - Dewars must arrive:
 - **Monday – Friday, except bank holidays**
 - **At least 1 day before beamtime**
 - **Before 12h00 (mid-day)**
 - Dewars must be properly labelled (especially those hand-delivered):
 - Beamline's Name
 - Local Contact's Name
 - Sender's Name & Telephone number
 - Please include:
 - Hard Disk Drive (≥ 1 TB)
 - **Return Label or Transporter's Return Air Way Bill**



Note: The PX2-A CATS Dewar holds nine (9) uni-pucks!

Preparing for REMOTE ACCESS: Account Activation

- Before your remote access session:
 - Obtain your SUNSET **project number and password**
 - E.g. 20190987, wR3nJH9s
 - Declare the remote participants
 - Declare the **Local Contact as the Session Responsible**
 - Test the connection to SOLEIL
 - Remote Access is activated **23 hours before** your session
 - and deactivates **25 hours after** your session
 - Contact the BL staff to set up a Remote Test
 - Send us a telephone number or Skype address




- IMPORTANT, you need:
 - To Install NoMachine
 - Also XQuartz for Mac
 - Decent internet connection
 - e.g. fiber optic, ADSL, 4G or 5G, not a modem
 - QWERTY keyboard
 - Try the command `setxkbmap fr` for AZERTY keyboards
 - Try the command `setxkbmap us` for QWERTY keyboards
 - Mouse with left & right buttons and scroll-wheel
 - To ALIGN, MOUNT and ROTATE samples
 - High resolution display(s) – **This is important!**
 - MXCuBE is optimised for **2560 x 1440**
 - **Otherwise, you will not see the bottom part of the MXCuBE window**



Testing your Internet Connection Speed

- Speedcheck.org
 - Latency
 - Return « ping » of small packets
 - Measure of Responsiveness
 - Download
 - Multiple connections opened
 - Send large file on all connections
 - Upload
 - Multiple connections opened
 - Receive large file on all connections



Connection Type	Latency [ms]	Download [Mbps]	Upload [Mbps]	Comments
Phone Hotspot @ PX2-A	120 - 250	0.03 - 0.43	0.33 – 1.00	5 s delays via Google Meet, sluggish but possible
SOLEIL Site @ PX2-A	50-51	120-129	127- 141	SOLEIL Site Network for office PCs, no connection to BL
SOLEIL Wifi @ PX2-A	22-39	144-178	168-192	Via SOLEIL Wifi Network on PX2-A (17/03/2023)
Freebox 4k @ chez Bill, Wifi < 30 cm	25-90	45-65	50-80	NX, MXCuBE, Albula, Webcams, GMeet, Zoom (on 2 nd Mac)...
Freebox 4k @ chez Bill, Wifi > 30 m*	35-45	15-55	10-70	Wifi station NOT in line of site, frequent cuts in connection


Updates and Changes

- SOLEIL has a new remote access server:
 - **Host** = `remote.synchrotron-soleil.fr`
 - **Port** = `4000`
 - **Protocol** = `NX`
 - `nx-vip.synchrotron-soleil.fr` should **NOT** be used
- PROXIMA 2A has a new MXCuBE workstation:
 - `proxima2a-pc4` (pc4)
 - Custom Display ONLY
 - `proxima2a-10` (p10) = old workstation
 - Still available via `remote.synchrotron-soleil.fr`
 - Physical and Custom Displays available
- NoMachine version 6
 - All workstations have been upgraded the NoMachine version 6
- ALBULA & ADXV
 - ADXV launches automatically with MXCuBE
 - A spot search algorithm (TIOGA) will display blue boxes around spots
 - ALBULA can be launched via command line on `proxima2a-pc4`
 - `albula` or `albula_4.1`


Connect via NoMachine (1): First Steps



- Start by clicking on the NoMachine icon
 - Within 23 hours before your session



- Click on remote.synchrotron-soleil.fr



- Enter your project-id and its password
 - For example: **20230123** [**w8d9kJ43**]



- Click on [proxima2a-pc4](#)
 - This is a new workstation that replaces **proxima2a-10**


Connect via NoMachine (2): Custom Display




- Click on: **Create a new desktop or custom session**




- Click on: **Create a new custom session**



- Click on **Continue**
 - Activate Run the console
 - Activate Run the command in a floating window






- Machine and Display options window
 - Do not modify the options
 - You will be connected soon...





- Success!
 - This last step may take awhile (> 60 seconds)

PX2-A Overview: Sample Environment & Performances



X-Terminals

- NoMachine should display:
 - A NoMachine window with Settings Icons, and
 - A purple gnome-terminal window with a prompt
 - If not, then a white X-terminal should display
 - Type: **gnome-terminal**
- If your key board is AZERTY
 - **setxkbmap fr** for AZERTY
 - **setxkbmap us** for QWERTY
- In the **purple** gnome-terminal , open 4 tabs:
 - **<shift-ctrl-t>**
 - Opens a new tab
 - **mxcube**
 - Opens an MXCuBE window
 - **albula**
 - Opens a window to visualise diffraction images
 - **firefox -P**
 - Opens a browser window for web cams
 - **ssh -X process1**
 - Connects to the data processing server
- To log out:
 - Type **exit** or **<ctrl-d>** in the X-terminals
 - This closes the session under NoMachine



MXCuBE Overview

- Type « **mxcube** » in a terminal window
 - Less than 2 minutes to start up
 - Logfile output is displayed in the terminal window
- Login to MXCuBE with your project number and password:
 - For example: **20230123** & **rF67wEr325**
 - Click **OK** when the small dialogue box appears

Goniometer positions & Action Buttons

Login

The screenshot displays the MXCuBE web interface with several key components highlighted by red dashed boxes and labels:


- OAV Microscope View:** A large central window showing a microscope image of a sample.
- Data Collection Methods:** A central panel containing various acquisition parameters such as 'Standard Collection', 'Acquisition', 'Exposure time (s)', 'Resolution (Å)', and 'Unit cell'.
- Sample Tree:** A panel on the right side showing a hierarchical list of sample positions (e.g., Puck 1, Puck 2) and their corresponding image numbers.
- Messages & Log:** A panel at the bottom center displaying system messages, including a notification that 'Data collection is enabled'.
- Beamline information:** A panel at the bottom right showing machine status, resolution (6.776 Å), current (800.00 mA), and energy (12.6500 keV).
- State information:** A panel at the bottom left showing the current state of the diffractometer and sample changer.

State information

Beamline information

Firefox & Webcams

- Open a tab and type:
 - **firefox** or **firefox -P**
 - Set AutoDetect Proxy settings
 - Preferences/Network Setting/ AutoDetect Proxy Settings
 - <http://cam14/view/view.shtml>
 - Four views of sample and beamline
 - <http://cam6/view/view.shtml>
 - View above hutch door
 - <http://cam8/view/view.shtml>
 - Overhead view of CATS Dewar
 - <http://cam13/view/view.shtml>
 - Motorized view of CATS, MD2 & EIGER
 - <http://cam1/view/view.shtml>
 - Motorized view of sample pin




Cam14 views

ALBULA


- Open a tab and type:

- **albula**


- Visualises diffraction images
- Click **Tools**
 - to change colour map, add resolution rings, etc...
- Click **Autoload/Eiger_Monitor**
 - To visualise the last image collected
 - Click on the DOWN ARROW to verify
 - » Host: 172.19.10.26
 - » Port: 80
 - » Min Pause (s): 1
- Click **File/Open**
 - Select a master file
 - » ***_master.h5**
 - Wedges of summed images
 - » ***_sum10_master.h5**
- Click on **Auto Contrast**
 - To adjust the contrast
- **Known Bugs**
 - Refresh the resolution rings with the scroll wheel
 - Differences between **pc4** and **p10**:
 - » AutoLoad does not work on **p10**
 - » Only SUM10 images on **p10**
 - » When remote, the **pc4** physical display needs to be logged in as **project_number**



ADXV




- **ADXV_FOLLOW**
 - Displays diffraction images automatically during collections
 - Launched with MXCuBE
 - If not, in an X-terminal window type:
 - **adxv_follow**
 - or **adxv** to display single images
 - TIOGA
 - Fast spot finding algorithm
 - Appears as blue squares
 - It may become “out-of-sync”
- **ADXV Commands**
 - Loads *.**cbf.gz** files
 - Video
 - Files = Sums N images
 - Stride = Skips N images
 - Type “**h**” to auto-adjust contrast




EIGER HDF5 File Architecture

- HDF5 files are “structured containers” made up of two types of files:
 - ONE **MASTER** file
 - MANY **DATA** files
- MASTER file
 - ***_master.h5**
 - Meta data
 - X-ray wavelength
 - Distance
 - Beamcenter
 - Goniometer angles
 - **Pointers to DATA files**
 - etc...
- DATA files
 - ***_data.h5**
 - Blocks of images per file
 - 1 to 1000s
- **DO NOT CHANGE HDF5 FILE NAMES!!!**
 - Keep master and data files in same directory
- Editing scripts
 - **fix_negative_transformation.py**
 - Fixes a geometrical parameter required by DIALS & XIA2



Data processing on “Process1”

- Open an X-terminal, type:
 - **ssh -X process1**
 - 288-core server
 - Dedicated to processing
- Some basic commands:
 - **going**
 - Go the latest data collection directory
 - **/nfs/data2/2023_Run2/20220789/2023-03-17/RAW_DATA/MyProj/process**
 - **goxdsme**
 - Go to the latest directory and launch XDSME on latest data collected
 - **/nfs/data2/2023_Run2/20220789/2023-03-17/RAW_DATA/MyProj/process/xdsme_MyXtal_1**
 - **goxdsme --brute**
 - Process for difficult indexing cases
 - **goxdsme --weak**
 - Process for weakly diffracting data
 - **xdsme -help**
 - More information...





Data processing with XDSME

- For **quick** processing in the X-terminal type:
 - **going**
 - To go to the « process » directory of the **most recently** collected data
 - **goxdsme**
 - To process the **most recently** collected data
- For processing older data :
 - Go to the directory of the **most recently** collected data
 - **going**
 - Change to the desired « process » directory, for example:
 - **cd ../../../../puck4356/pin15/process**
 - Process the specific master file with **xdsme** (not goxdsme):
 - **xdsme ../pin15_2_master.h5**
- Common XDSME options:
 - **xdsme --brute ../pin15_2_master.h5**
 - Tries hard to index...
 - **xdsme --weak ../pin15_2_master.h5**
 - For cases with weak diffraction
 - **xdsme -r 3.5 ../pin15_2_master.h5**
 - Limits the resolution to 3.5 Å
 - **xdsme -s P2 -c "102 124 87 90 97 90" ../pin15_2_master.h5**
 - Inputs space group and unit cell
 - **xdsme -F 301 -L 3300 ../pin15_2_master.h5**
 - Processes from images 301 (first) to 3300 (last)
 - **xdsme -h**
 - Prints out "help" and more options

Help, Hints, FAQs & Bugs

- Some DON'Ts
 - DON'T use VPN at the same time as No Machine
 - Disconnect your VPN software (Junos Pulse, etc...)
- Some DOs
 - DO use GNOME-TERMINAL
 - Use the command: **gnome-terminal**
 - DO create your own FIREFOX profile
 - Use the "-P" option: **firefox -P**
- Keyboard Settings
 - Use the command **setxkbmap fr** for AZERTY
 - Use the command **setxkbmap us** for QWERTY
- NoMachine Errors
 - Ooops! Negotiation failed...
 - Your user_id does not have permission to connect to the NX server
 - DO NOT USE: **nx-vip.synchrotron-soleil.fr**
 - INSTEAD USE: **remote.synchrotron-soleil.fr**
 - Ooops! The remote host "name" can not be found...
 - Check the host name, port number and protocol are correct
 - Authentication failed...
 - Check the project_id and password are correct
 - Ooops! Could not connect to server. Error 60 Timed out
 - The project_id may not have permission to connect
 - Verify the current time is within 23 hours before and 25 hours after the session time
- HINTS
 - If the Data Collection widgets are GREYED OUT
 - Check that you are logged in to MXCuBE
 - Click on the Sample Line in the Sample Tree
 - To highlight which sample is mounted
 - Use the "ISpyB" button to refresh the Sample Tree
 - If you can NOT find your HDF5 images
 - Check the project_number in the pathname
 - Check that "downloader" is running under your user_id
 - Type in an X-terminal: **downloader status**
 - If it is NOT RUNNING : **downloader start**
 - If it is RUNNING under another user_id
 - » Call the Local Contact
 - Robot collisions
 - Check the position of the robot arm via the webcams
 - Try the **SAFE** button in MXCuBE
 - If not, the call the BL or EHOs (ext 9797)
 - If the microscope image is dark (deactivated)
 - Check that the "camera" device server is running
 - Type in an X-terminal: **camera restart**
- Telephone numbers
 - Beamline telephone numbers:
 - +33-1-6935-8181 (wireless)
 - +33-1-6935-8180 (landline)
 - Experimental Hall Operators (EHOs)
 - +33-1-6935-9797 (after 11 pm)



NoMachine Connection Errors



“Oops! The remote host can not be found”:

- 1) Click EDIT
- 2) Check the host, port and protocol
 - Look for typographical errors

Host = remote.synchrotron-soleil.fr
Port = 4000
Protocol = NX



- 1) The username/password is wrong, or
- 2) REMOTE ACCESS has **not** been activated for the project number.
 - Please wait until 23 hours before the session
 - If not then contact the LC




“Oops! Could not connect to the server”

- 1) Try switching off any VPN?
- 2) Maybe the NX node is saturated?

- Please contact the LC

?



“Oops! The session negotiation failed.”

- 1) Maybe too early to connect?
- 2) Maybe the NX node is saturated?

- Please contact the LC

?

