

# LUCIA

beamline

**XRF, XANES and EXAFS**

ORGANISMS



TISSUES



CELLS



ORGANELLES



COMPLEXES



PROTEINS




ATOMS



# SAMPLES

## TYPES

- Cryogenic or freeze-dried thin sections from vegetal, animal, fungi's tissues. Thickness between 5 and 50  $\mu\text{m}$ .
- Cryo-fixed or freeze-dried cells.
- Freeze-dried grounded material from living organism or from soils, sediments, suspended matters,... pressed as pellets.



## SAMPLE HOLDERS

- Thin section placed between two Ultralene foils in Cu-sample holder.
- Cells cultivated onto Ultralene or  $\text{Si}_3\text{N}_4$  windows.
- Pellets placed on Cu-sample holder.



## ENVIRONMENT

- **Vacuum** chamber.
- **Cryogenic conditions** using  $\text{N}_2$  or He cryostat.
- No glass or quartz slide. Resin impregnation is not recommended.
- Constraints on the sample holder and sample preparation may vary as a function of the selected energy, the nature and concentration of the elements of interest and of the matrix, and of the type of analysis.

# TECHNIQUE

LUCIA is a **microfocused beamline** in the 0.8-8 keV energy range. The spatial stability of the beam spot over a wide spectral range enables elemental distribution studies by **micro-fluorescence X-ray**



**spectroscopy** ( $\mu$ -XRF) as well as elemental speciation by **X-ray absorption spectroscopy** (XANES and EXAFS) to determine the speciation (local chemistry, quantitative determination of the local geometric structure around the absorbing atom) in heterogeneous samples.

## Micro-fluorescence X-ray spectroscopy ( $\mu$ XRF)



High spatial resolution mapping ( $2.5 \mu\text{m}$ ) on large fields of view (up to  $20 \times 20 \text{ mm}^2$ ) with excitation energies between 0.8 and 8 keV.



Fast scans (FlyScan) to analyze a large number of samples.

## X-ray absorption spectroscopy



XANES (X-ray Absorption Near-Edge structure) and EXAFS (Extended X-ray absorption fine structure) at the K-edge of elements from Na to Fe, L edges from Ni to Gd, and M edges of rare earths and actinides.

## Combination of micro-XRF and micro-XANES




After a  $\mu$ XRF cartography of the sample, interesting spots can be analyzed by  $\mu$ XAS to determine locally the speciation of the elements and how this depends on the different elemental composition.



# HIGHLIGHTS

$\mu$ XRF mapping was used to determine the tissue-level distribution of Mn, Ca, and other elements in different tissues of *G. Meisner*, a Mn-accumulator plant. Mn concentrations were highest in leaf epidermal tissues, in cortex and vascular tissues of primary roots. The co-localization of Ca and Mn might in some way facilitate the formation of mixed Ca-Mn oxides

C. Bihanic et al., **Sci. Rep.** (2021)



Combination of  $\mu$ XRF and  $\mu$ XANES at **Ag L<sub>3</sub>-edge** evidenced the physico-chemical changes of Ag colloidal suspension following ingestion and digestion by mice.



The XRF map of a cryo-thin section of mouse colonic content displays a co-localization of Ag and S.  $\mu$ XANES evidence the formation of Ag complexes with thiol groups of proteins and/or peptides.

K. Gillois et al., **Environ. Sci. Nano.** (2023)



# REFERENCES

- D. Vantelon et al., **The LUCIA beamline at SOLEIL**. J. Synchrotron Rad., (2016)
- A.M. Flank et al., **LUCIA, a microfocus soft XAS beamline**. Nucl Instrum Methods Phys Res B, (2006)



**More information on LUCIA publications web page**

## COMPLEMENTARY BEAMLINES

**NANOSCOPIUM** for quantitative imaging (morphology, elemental composition and chemical speciation) from nano to microscopic size. **ANATOMIX** to obtain two- and three dimensional radiographic images of bulk volume samples of macroscopic and nano size.

**SMIS**, to obtain images using infrared spectromicroscopy. **DISCO** to probe the fluorescence of living and fixed samples using VUV to visible light.



**LUCIA end-station**

# CONTACT

**Delphine Vantelon**

*Beamline Manager*

✉ delphine.vantelon@synchrotron-soleil.fr

☎ +33 (0)1 69 35 96 94

**Camille Rivard**

*Beamline scientist for biology projects*

✉ camille.rivard@synchrotron-soleil.fr

☎ +33 (0)1 69 35 94 71

## Health & Well-Being at SOLEIL



Link to the  
web page

🐦 @HelioBio\_SOILEIL

SOLEIL's Health and Well-being Scientific Section is composed of 30 scientific experts from different fields. Through collaborative and science-driven approaches, the Section offers the community a coherent portfolio of state-of-the-art techniques to serve scientific and societal health-related challenges.



L'Orme des Merisiers - Départementale 128  
91190 Saint-Aubin - FRANCE  
[www.synchrotron-soleil.fr](http://www.synchrotron-soleil.fr)

