

# SMIS beamline

**MICRO**

ORGANISMS



TISSUES



CELLS



ORGANELLES



COMPLEXES



PROTEINS



ATOMS



# SAMPLES

## TYPES

- **Cells (fixed):** 2D cultures on infrared transparent windows (CaF<sub>2</sub>, ZnS...), primary cells, cell lines, yeast, bacteria, protozoa...
- **Tissues:** thin sections (2-10 μm thick , 20-30 μm for vegetal), 3D cell cultures, cryo-sectioned without embedding, deposited on IR transparent windows
- **Cells (living ):** require a microfluidic device with temperature control (available at the beamline)

## ENVIRONMENTS

- **Transparent substrates:** compatible for cell growth, and/or with different spectral domains
- **Temperature stage:** -180 to 600°C, purged (Linkam FTIR600)
- **Microfluidic device:** CaF<sub>2</sub> windows, 2 intake and 2 exit ports, heating elements, temperature control
- **Stretching device:** with temperature control (Linkam TFT350), 2 and 20 N force transducers

## OFFLINE INSTRUMENTS

- **Raman microscope:**  
TFS DXR, 532, 633 and 780 nm lasers, 50x and 100x objectives, 1 μm resolution
- **Imaging IR microscope:**  
Cary 620, 128x128 pixel Focal Plane Array detector to image large samples (mm<sup>2</sup> to cm<sup>2</sup> size) quickly with 20-30 μm resolution.



# TECHNIQUE

**SMIS** is an infrared microspectroscopy beamline that is organized in two branches: **SMIS Micro** and **SMIS Nano**. The SMIS Micro branch specializes in synchrotron infrared microspectroscopy and imaging for life science, cultural heritage, material under extreme conditions, polymers, astrophysics...


The **SMIS micro** branch operates two confocal infrared microscopes coupled to synchrotron to probe biological samples at the cellular and subcellular resolution (3-15  $\mu\text{m}$ )




Infrared microspectroscopy can be used to:

- ① **Identify biomarkers** such as biomolecules or minerals within cells and tissues.
- ② **Quantify** the concentration of known biomarkers (polysaccharides, lipids, proteins, metabolites...).
- ③ Estimate the **conformation** of proteins in cells and tissues, detect amyloids plaque and inclusions; lipid conformation...
- ④ **Fingerprinting** cells and tissues to identify and delineate **cell types**, **metabolic** or **pathologic** states.

**Identification and quantification of creatin** in rat brain hippocampus




**Conformation of proteins in inclusions in Huntington Disease patient brain**




# HIGHLIGHTS

Chronic Myeloid Leukemia (CML) can be treated efficiently by Tyrosine Kinase Inhibitors (TKI). Resistance to TKI therapies may develop during treatment due to mutations occurring in the ABL-kinase domain of the BCR-ABL oncogene. This has led to the development of a second and a third generation of TKI drugs. The T315I mutation in the bcr-abl gene renders leukemic cells resistant to all three generations of TKIs.

Synchrotron  $\mu$ FTIR can be used to measure the chemical composition of individual myeloid progenitor cells and extract a spectral signature of single cells expressing the T315I mutation.



Differentiation of murine embryonic stem cells expressing WT or T315I bcr-abl by SR- $\mu$ FTIR and PCA. Peaks at 1242, 1122, 1088, 966 cm<sup>-1</sup> are tentatively associated to RNA ribose-phosphate backbone.



Extinction of the T315I signature in the doxinducible UT7IND cell lines carrying WT or T315I bcr-abl.

# REFERENCES

- ② P. Dumas, et al. Synchrotron infrared microscopy at the French synchrotron facility SOLEIL Infrared Phys. Technol. 49(1-2), 152 (2006).  
10.1016/j.infrared.2006.01.030
- ② P. Dumas, et al., Adding synchrotron radiation to infrared microspectroscopy: what's new in biomedical applications?, Trends Biotechnol. 25(1), 40 (2007).  
10.1016/j.tibtech.2006.11.002.



**More information on SMIS publications web page**

## COMPLEMENTARY BEAMLINES

**ANATOMIX:** obtain three-dimensional X-ray tomography images of bulk volume samples at microscopic resolution.

**NANOSCOPIUM:** micro to nano morphology, elemental composition and chemical speciation.

**DISCO:** chemical imaging using auto-fluorescence microspectroscopy.

**LUCIA:** X-ray microprobe ( $\mu$ -XAS,  $\mu$ -XRF).



# CONTACT

**Ferenc Borondics**

*Beamline Manager*

✉ ferenc.borondics@synchrotron-soleil.fr

☎ +33 (0)1 69 35 81 92

**Christophe Sandt**

*Beamline scientist, life sciences*

✉ christophe.sandt@synchrotron-soleil.fr

☎ +33 (0)1 69 35 81 07

## Health & Well-Being at SOLEIL



Link to the  
web page

SOLEIL's Health and Well-being Scientific Section is composed of 30 scientific experts from different fields. Through collaborative and science-driven approaches, the Section offers the community a coherent portfolio of state-of-the-art techniques to serve scientific and societal health-related challenges.



L'Orme des Merisiers - Départementale 128  
91190 Saint-Aubin - FRANCE  
[www.synchrotron-soleil.fr](http://www.synchrotron-soleil.fr)